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When it is desirable to exhaust air to the atmosphere, it should be via a 100% exhaust system (i.e., a system that does not recirculate its exhaust air into other parts of the building). The recommended exhaust system connection for types A1 and A2 cabinets is an exhaust canopy connection as is shown in annex E, figures E2 and E3. Every canopy design must be tested to determine the airflow rate exhausted by the canopy that will assure performance. Whenever the cabinet is field certified, the minimum exhaust flow by the canopy should be verified by measurements using the approved instruments and techniques cited in annex A, sections A.9 and A. 10. No type A cabinet should ever be hard connected to an exhaust system (see figures E4 and E5).

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When Type A1 and A2 cabinets are found to be directly attached to the exhaust system and vented to the outside without the use of an exhaust canopy, it is recommended that the exhaust connection be modified to an exhaust canopy.

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Types B1 and B2 cabinets are to be vented outside the building without recirculation. The venting system should include a leak-tight duct a damper in the duct near the cabinet to permit flow adjustment closure and decontamination, and an external exhaust fan as the final system component (see annex E, figure E6). The exhaust fan should be sized to deliver the required exhaust airflow (as specified by the cabinet manufacturer), considering pressure losses in the duct and allowing at least 2 in. w. g. (500Pa) for a dirty HEPA filter. If a charcoal filter is used downstream of the HEPA filter, an additional pressure capacity equal to the manufacturer's recommended resistance should be provided. An alarm should be provided at the cabinet to indicate loss of exhaust flow. This can be a exhaust volume flow measuring device in the duct downstream of the exhaust filter, sail switch at the fan discharge, or flow measuring station in the exhaust duct. It is recommended that each Type B1 or B2 cabinet have its own (dedicated) exhaust system. The cabinet should be interlocked with the blower in the duct or the building system to prevent pressurization of the exhaust system. In addition, cabinets hard connected to an exhaust system should not be turned off.

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